

IonSense CERAMIC TUBE LENGTHS FOR THE VAPUR INTERFACE DART 100, DART-ET, 1D DART-SVP, "3+D" DART-SVP

| Manufacturer | Flange Version | LC/MS System | DART 100 | DART-ET * | 1D DART-SVP | 2D DART-SVP |
|------------------------------|-----------------------|---------------------|----------------------|-----------------|----------------|-----------------|
| Agilent ID 3.18 mm | JAN-10 - Current | All Agilent Systems | Length 78 mm | Length 98 mm | Length 83 mm | Length 167 mm |
| | | | | | | |
| Bruker ID 4.75 mm | JUL-09 - Current | Bruker | Length 37 mm | Length 65 mm | Length 50 mm | Length 134 mm |
| | | | | | | |
| ABI/MDS ID 4.75 mm | Prior to DEC-09 | SCIEX API 100 LC | Length 127 mm (5 in) | Length 147 mm | Length 132 mm | Length 216 mm |
| | Current | SCIEX 3000 | | | | |
| | Current | QTRAP | Length 98 mm | Length 114.5 mm | Length 99.5 mm | Length 183.5 mm |
| | DEC-09 - Current | SCIEX 4000 | Length 58 mm | Length 79 mm | Length 64 mm | Length 148 mm |
| | DEC-09 - Current | SCIEX 5000 | Length 58 mm | Length 79 mm | Length 64 mm | Length 148 mm |
| JEOL ID 4.75 mm | Prior to JAN-10 | AccuTOF | Length 58 mm | Length 78 mm | Length 63 mm | Length 147 mm |
| | JAN-10 - Current | AccuTOF | Length 35 mm | Length 56 mm | Length 41 mm | Length 125 mm |
| | | | | | | |
| Thermo ID 4.75 mm | Prior to JAN-10 | Classic/DECA | Length 68 mm | Length 88 mm | Length 73 mm | Length 157 mm |
| | JAN-10 - Current | Classic/DECA | Length 73 mm | Length 94 mm | Length 79 mm | Length 163 mm |
| | Current | DECA XP | | | | |
| | Prior to DEC-09 | LTQ & Exactive | Length 42 mm | Length 68 mm | Length 53 mm | Length 137 mm |
| | DEC-09 - Current | LTQ & Exactive | | Length 57 mm | Length 42 mm | Length 126 mm |
| Current | TSQ 7000 | | | | | |
| Waters ID 4.75 mm | Current | Fishbowl, ZQ, TOF | Length 42 mm | Length 62 mm | Length 47 mm | Length 131 mm |
| | Current | Quattro | Length 42 mm | Length 62 mm | Length 47 mm | Length 131 mm |
| | Current | Premiere | Length 42 mm | Length 62 mm | Length 47 mm | Length 131 mm |
| | Coming in 2010 | Xevo | | | | |

* w/ SVP Exp. Blocks (NOV-09 - Pres.)